

## Access Array™ Target-specific Primers

Kit Part Numbers: ASY-AA, ASY-AAX

### Intended Use

Access Array Target-specific primers are intended for use with the Access Array System, Barcode Primer Plates, and Sequencing Primer Kit for next-generation sequencing experiments. Please refer to the *Access Array System for 454 Sequencing Platform User Guide*, PN 68000158, and *Access Array System for Illumina Sequencing Platform User Guide*, PN 100-3770, for protocols.

Contents				
PART NUMBER	REAGENT	CONCENTRATION of each individual primer*	VOLUME	STORAGE
ASY-AA	Mixed	50 µM	100 µL	-20°C
ASY-AAX	Forward and Reverse Primers	60 µM	100 µL	

\*ASY-AA-WET plate contains 50 µL of 100 µM forward and 50 µL of 100 µM reverse primers

\*ASY-AAX plate contains 50 µL of 120 µM forward and 50 µL of 120 µM reverse primers.

### Handling:

Before using reagents:

1. Warm to room temperature.
2. Centrifuge for 30 seconds to collect materials.
3. Remove original plate seal. Original seal is intended for single use only.
4. Manually mix the contents of the plate. (Do not vortex to avoid cross contamination.)
5. Reseal the plate with a new seal.
6. Centrifuge for 30 seconds to collect contents.

### Safety:

Use standard laboratory safety protocols. For MSDS information, contact Fluidigm Technical Support.

### Technical Support:

United States: 1-866-FLUIDLINE (1.866.358.4354) or techsupport@fluidigm.com

Outside the United States: 650-266-6100 | [www.fluidigm.com/support.htm](http://www.fluidigm.com/support.htm)

Fluidigm, the Fluidigm logo, and Access Array are trademarks or registered trademarks of Fluidigm Corporation in the U.S. and/or other countries. Fluidigm recommends you only use licensed PCR assay reagents from authorized sources. ©Fluidigm Corporation. All rights reserved. For research use only. Not for diagnostic use.