

EQ[™] Four Element Calibration Beads

WARNING Before handling any chemicals, refer to the Safety Data Sheet (SDS) provided by the manufacturer, and observe all relevant precautions.

Catalog#: 201078

Contents and Storage:

- EQ[™] Four Element Calibration Beads contain natural abundance cerium (140/142Ce), europium (151/153Eu), holmium (165Ho), and lutetium (175/176Lu).
- 100 mL of beads at 3.3x10E+05 beads/mL are provided.
- Store at 4 °C (stable for at least 6 months).

Directions for Use

For data normalization samples are re-suspended in 0.1X EQ Beads and acquired according to the following procedure:

- 1 Shake beads vigorously** to re-suspend.
- 2 Prepare sufficient volume of 0.1X beads to re-suspend all samples in the experiment by diluting 1 part beads to 9 parts deionized water (18MΩ, ultrapure).
- 3 After completing the final water wash of the cell sample, but immediately prior to injecting sample for analysis, re-suspend the resultant cell pellet in 0.1X bead solution (~33,000 EQ beads per mL).
- 4 Filter the cell-bead suspension through a 35 to 45 μm mesh immediately prior to injection into the instrument.
- 5 When creating the instrument Template for data collection, select Ce140, Eu151, Eu153, Ho165 and Lu175 as channels to be collected (in addition to other channels necessary for the experiment).
- 6 Follow the directions in the User Guide, "Normalization of Mass Cytometry Data using EQ[™] Four Element Beads" for additional information on collection and analysis of data for normalization.

For technical support visit fluidigm.com/support

EMAIL United States and other countries not in Europe or Asia support.northamerica@fluidigm.com | Europe support.europe@fluidigm.com | Asia techsupportasia@fluidigm.com
PHONE United States (toll-free) +1 866 358 4354 | Europe +33 1 60 92 42 40 | Japan +81 3 3662 2150 | China +86 21 3255 8368 | All other countries +1 650 266 6100

For Research Use Only. Not for use in diagnostic procedures.

Information in this publication is subject to change without notice. Use standard laboratory safety protocols. **Safety data sheet information** fluidigm.com/sds
Patent and license information fluidigm.com/legalnotices **EU's WEEE directive information** fluidigm.com/compliance | Fluidigm, the Fluidigm logo, CyTOF, EQ, and Maxpar are trademarks or registered trademarks of Fluidigm Corporation in the United States and/or other countries. All other trademarks are the sole property of their respective owners. © 2015 Fluidigm Corporation. All rights reserved. 04/2015