

Advanta Immuno-Oncology Gene Expression Assay Panel B Assay Plate

Kit Part Number: 101-6146

Purpose

This kit provides necessary assays for immuno-oncology gene expression with 96.96 Dynamic Array™ IFCs (integrated fluidic circuits) for gene expression. The kit contains 79 assays in one 96-well plate.

Contents

Part Number	Component	Stock Concentration	Volume	Storage
101-5851	Advanta™ Immuno-Oncology Expression Assay Panel B assay plate	20X	10 µL/assay	-20 °C

Handling

Before using assay, refer to the Advanta Immuno-Oncology Gene Expression Assay (PN 101-6108).

Safety



DANGER Advanta Immuno-Oncology Gene Expression Assay Panel B assay plate contains formamide. May damage fertility or the unborn child if swallowed. Suspected of causing cancer if swallowed. If exposed or concerned: Get medical advice/attention. Do not handle until all safety precautions have been read and understood. Obtain special instructions before use. Wear protective gloves/protective clothing/eye protection/face protection. Read the safety data sheet.

Use standard laboratory safety protocols. Read and understand the safety data sheets (SDSs) before handling chemicals. To obtain SDSs, go to fluidigm.com/sds and search for the SDS using either the product name or the part number.

For technical support visit fluidigm.com/support.

North America +1 650 266 6100 | Toll-free (US/CAN): 866 358 4354 | techsupport@fluidigm.com **Latin America** +1 650 266 6100 | techsupportlatam@fluidigm.com
Europe/Middle East/Africa/Russia +44 1223 859941 | techsupporteurope@fluidigm.com **China (excluding Hong Kong)** +86 21 3255 8368 | techsupportchina@fluidigm.com
Japan +81 3 3662 2150 | techsupportjapan@fluidigm.com **All other Asian countries/India/Australia** +1 650 266 6100 | techsupportasia@fluidigm.com

For Research Use Only. Not for use in diagnostic procedures.

NOTICE TO PURCHASER: LIMITED LICENSE: Use of this product is covered by one or more US patents. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim, no right to perform any patented method, and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting outlicensing@lifetech.com or Out Licensing, Thermo Fisher Scientific, 5781 Van Allen Way, Carlsbad, California 92008.

Information in this publication is subject to change without notice. Patent and license information: fluidigm.com/legalnotices. Fluidigm, the Fluidigm logo, Advanta, and Dynamic Array are trademarks or registered trademarks of Fluidigm Corporation in the United States and/or other countries. © 2017 Fluidigm Corporation. All rights reserved. 09/2017