

Announcing New Maxpar Reagents for Mass Cytometry

We're pleased to announce the addition of 62 new metal-conjugated antibodies to the Fluidigm Maxpar[®] reagent catalog: 19 against mouse epitopes, 39 against human epitopes and 4 that are cross-reactive.

Scroll down to view the complete list. To place an order, [contact your sales representative](#).

May the discoveries continue.

Sincerely,

The Customer Support Team

Maxpar Reagent Catalog

NOTE Catalog numbers ending in A denote a 50-test size, while those ending in B denote a 100-test size. All reagents can be ordered in a 25-test size by replacing the last letter in the catalog number with the letter C.

Anti-Human	Cat. No.	Application
CD44s Pan Specific (691534)— ¹⁵³ Eu	3153021B	CD and Surface Markers
CD49a (TS2/7)— ¹⁶³ Dy	3163015B	CD and Surface Markers
CD49c (ASC-1)— ¹⁶¹ Dy	3161016B	CD and Surface Markers
CD54 (HA58)— ¹⁷⁰ Er	3170014B	CD and Surface Markers
CD57 (HCD57)— ¹⁷⁶ Yb	3176019B	CD and Surface Markers
CD94 (HP-3D9)— ¹⁷⁴ Yb	3174015B	CD and Surface Markers
CD100 (A8)— ¹⁵⁶ Gd	3156021B	CD and Surface Markers
CD127 (A019D5)— ¹⁶⁸ Er	3168017B	CD and Surface Markers
CD130 (2E1B02)— ¹⁶⁰ Gd	3160017B	CD and Surface Markers
CD134/OX40 (ACT35)— ¹⁵⁰ Nd	3150023B	CD and Surface Markers
CD181 (8F1/CXCR1)— ¹⁴² Nd	3142009B	CD and Surface Markers
CD182 (5E8/CXCR2)— ¹⁴⁷ Sm	3147010B	CD and Surface Markers

Anti-Human	Cat. No.	Application
CD193/CCR3 (5E8)— ¹⁷⁵ Lu	3175025B	CD and Surface Markers
CD223/LAG3 (874501)— ¹⁶⁵ Ho	3165028B	CD and Surface Markers
CD235 (HIR2)— ¹⁷⁵ Lu	3175027B	CD and Surface Markers
CD257 (T7-241)— ¹⁶³ Dy	3163013B	CD and Surface Markers
CD271 (C40-1457)— ¹⁴⁹ Sm	3149017B	CD and Surface Markers
CD274/PD-L1 (29E.2A3)— ¹⁵⁶ Gd	3156026B	CD and Surface Markers
CD278/ICOS (DX29)— ¹⁴³ Nd	3143020B	CD and Surface Markers
CD279/PD-1 (EH12.2H7)— ¹⁵⁵ Gd	3155009B	CD and Surface Markers
CD301/CLEC10a (H037G3)— ¹⁶⁷ Er	3167015B	CD and Surface Markers
CD325/N-Cadherin (8C11)— ¹⁴³ Nd	3143016B	CD and Surface Markers
B2-microglobulin (2M2)— ¹⁷⁰ Er	3170016B	CD and Surface Markers
DNMT3B (832121)— ¹⁶⁴ Dy	3164021B	Signaling and Transcription
Galanin (581403)— ¹⁴⁸ Nd	3148016B	Signaling and Transcription
IL-13 (JES10-5A2)— ¹⁶⁹ Tm	3169016B	Cytokines
IL-2 (MQ1-17H12)— ¹⁴⁴ Nd	3144021B	Cytokines
Integrin $\alpha 9\beta 1$ (Y9A2)— ¹⁶⁸ Er	3168013B	CD and Surface Markers
Invariant NK T (6B11)— ¹⁷⁰ Er	3170015B	Signaling and Transcription
LCK (LCK-01)— ¹⁵³ Eu	3153017B	Signaling and Transcription
Nestin (196908)— ¹⁴⁶ Nd	3146015B	Structure and Cytoskeleton
p21 Waf1/Cip1 (12D1)— ¹⁵⁹ Tb	3159026A	Cell Cycle and Proliferation
Progesterone Receptor (D8Q2J)— ¹⁴⁵ Nd	3145011A	Signaling and Transcription
Semaphorin 6A (169203)— ¹⁶³ Dy	3163018B	Signaling and Transcription
Siglec 8 (7C9)— ¹⁶⁴ Dy	3164017B	CD and Surface Markers
SR-AI/MSR (351615)— ¹⁵⁶ Gd	3156025B	Signaling and Transcription
Stat3 (124H6)— ¹⁷³ Yb	3173003A	Signaling and Transcription
TIGIT (MBSA43)— ¹⁵⁴ Sm	3154016B	CD and Surface Markers

Anti-Human	Cat. No.	Application
------------	----------	-------------

TrkA (165131)— ¹⁴⁹ Sm	3149019B	Signaling and Transcription
----------------------------------	----------	-----------------------------

Anti-Mouse	Cat. No.	Application
------------	----------	-------------

BATF (D7C5)— ¹⁵⁴ Sm	3154012A	Signaling and Transcription
--------------------------------	----------	-----------------------------

CD45RB (C363.16A)— ¹⁴⁵ Nd	3145012B	CD and Surface Markers
--------------------------------------	----------	------------------------

CD49d (R1-2)— ¹⁵¹ Eu	3151016B	CD and Surface Markers
---------------------------------	----------	------------------------

CD5 (53-7.3)— ¹⁴⁶ Nd	3146012B	CD and Surface Markers
---------------------------------	----------	------------------------

CD54/ICAM-I (YN1/1.7.4)— ¹⁶³ Dy	3163020B	CD and Surface Markers
--	----------	------------------------

CD66a/CEACAM1 (CC1)— ¹⁵³ Eu	3153022B	CD and Surface Markers
--	----------	------------------------

CD154/CD40 Ligand (MR1)— ¹⁷⁰ Er	3170011B	CD and Surface Markers
--	----------	------------------------

CD180/RP105 (RP/14)— ¹⁶² Dy	3162031B	CD and Surface Markers
--	----------	------------------------

CD185/CXCR5 (614641)— ¹⁴² Nd	3142008B	CD and Surface Markers
---	----------	------------------------

CD223/LAG3 (C9B73)— ¹⁷⁴ Yb	3174019B	CD and Surface Markers
---------------------------------------	----------	------------------------

CD272/BTLA (6F7)— ¹⁵⁶ Gd	3156028B	CD and Surface Markers
-------------------------------------	----------	------------------------

CD278/ICOS (7E.17G9)— ¹⁷⁶ Yb	3176014B	CD and Surface Markers
---	----------	------------------------

CD279/PD1 (29F.1A12)— ¹⁵⁹ Tb	3159024B	CD and Surface Markers
---	----------	------------------------

CD279/PD1 (J43)— ¹⁵⁹ Tb	3159023B	CD and Surface Markers
------------------------------------	----------	------------------------

CD357/GITR (DTA1)— ¹⁴³ Nd	3143019B	CD and Surface Markers
--------------------------------------	----------	------------------------

CD366/TIM3 (RMT3-23)— ¹⁶² Dy	3162029B	CD and Surface Markers
---	----------	------------------------

iNOS (CXNFT)— ¹⁶¹ Dy	3161011B	Signaling and Transcription
---------------------------------	----------	-----------------------------

RORgt (B2D)— ¹⁵⁹ Tb	3159019B	Signaling and Transcription
--------------------------------	----------	-----------------------------

TNFalpha (MP6-XT22)— ¹⁴¹ Pr	3141013B	Cytokines
--	----------	-----------

Cross-Reactive	Cat. No.	Application
Anti-Human/Mouse BCL-6 (K112-91)— ¹⁶³ Dy	3163012B	Signaling and Transcription
Anti-Mouse/Rat CD40 (HM40-3)— ¹⁶¹ Dy	3161020B	CD and Surface Markers
Pan-Actin (D18C11)— ¹⁷⁵ Lu	3175026A	Structure and Cytoskeleton
Vimentin (RV202)— ¹⁵⁶ Gd	3156023A	Structure and Cytoskeleton

For Research Use Only. Not for use in diagnostic procedures.

For technical support, visit fluidigm.com/support.

CORPORATE HEADQUARTERS
 7000 Shoreline Court, Suite 100
 South San Francisco, CA 94080 USA
 Toll-free: +1 866 359 4354 | Fax: +1 650 871 7152
fluidigm.com

SALES
 North America | +1 650 266 6170 | info-us@fluidigm.com
 Europe/EMEA | +33 1 60 92 42 40 | info-europe@fluidigm.com
 China (excluding Hong Kong) | +86 21 3255 8368 | info-china@fluidigm.com
 Japan | +81 3 3662 2150 | info-japan@fluidigm.com
 All other Asian countries | +1 650 266 6000 | info-asia@fluidigm.com
 Latin America | +1 650 266 6000 | info-latinamerica@fluidigm.com

Information in this publication is subject to change without notice. **Patent and license information:** fluidigm.com/legalnotices.
 Fluidigm, the Fluidigm logo and Maxpar are trademarks or registered trademarks of Fluidigm Corporation in the United States and/or other countries. © 2015 Fluidigm Corporation. All rights reserved. 09/2015