

Callisto System Software

v3.8.13 Enhancements

- **Networked use of the Callisto™ system.** New functionality to communicate with Callisto to modify experiments.
- **Improved run control.** New run control screen to modify experiments while they are running.
- **In-progress experiment modifications.** Modifying an experiment plan no longer requires stopping the experiment.
- **Dispensing map display.** New dispensing map and list of required reagents are displayed before running the experiment.
- **Calendar and schedule display.** New Calendar and Plan screens display the required attendance.
- **Step skipping and repeat.** Steps can be skipped or repeated on the Steps screen.

v3.8.13 Bug Fixes

- **Minimized software crashes.** Reduced memory allocation to prevent out of memory crashes.

For technical support visit fluidigm.com/support.

North America +1 650 266 6100 | Toll-free: 866 358 4354 in the US | techsupport@fluidigm.com Europe +33 1 60 92 42 40 | techsupporteurope@fluidigm.com

China (excluding Hong Kong) +86 21 3255 8368 | techsupportchina@fluidigm.com Japan +81 3 3662 2150 | techsupportjapan@fluidigm.com

All other Asian countries +1 650 266 6100 | techsupportasia@fluidigm.com Central and South America +1 650 266 6100 | techsupportlatam@fluidigm.com

For Research Use Only. Not for use in diagnostic procedures.

Information in this publication is subject to change without notice. Safety data sheet information: fluidigm.com/sds. Patent and license information: fluidigm.com/legalnotices. EU's WEEE directive information: fluidigm.com/compliance. | Fluidigm, the Fluidigm logo, and Callisto are trademarks or registered trademarks of Fluidigm Corporation in the United States and/or other countries. © 2016 Fluidigm Corporation. All rights reserved. 09/2016