

## Load Coil

Catalog#: 101803  
Package Size: 1 unit  
Instrument: CyTOF

## Technical Information

### Description:

The CyTOF® Load Coil is a copper wire coiled around the Torch Body and generates a radio frequency-induced electromagnetic field responsible for igniting and sustaining the Argon plasma. The correct coil pitch is maintained by the black coil core which is removed following installation.



## Contact Information:

Sales: [sales@DVSciences.com](mailto:sales@DVSciences.com) | Support: [Support.NorthAmerica@DVSciences.com](mailto:Support.NorthAmerica@DVSciences.com) or [Support.Europe@DVSciences.com](mailto:Support.Europe@DVSciences.com)  
[www.DVSciences.com](http://www.DVSciences.com) | For assistance by phone: +1-855-DVS-CYTO

**This product is intended for Research Use Only. This product is not intended for therapeutic or diagnostic purposes in humans or animals.**

**Vendor info & description quote as needed.**

List supplier of raw material antibody here only if targets/clones are unique, licensed, and/or indicated as required in our agreement.

© 2011 DVS Sciences Inc. The trademarks mentioned herein are the property of DVS Sciences Inc. or their respective owners.

DVS Sciences Inc. products are warranted to meet stated product specifications and to conform to label descriptions when used and stored properly. Unless otherwise stated, this warranty is limited to 6 months from date of sales for products used, handled and stored according to DVS Sciences, Inc. instructions. DVS Sciences, Inc. sole liability for the product is limited to replacement of the product or refund of the purchase price. Product may not be resold, modified for resale or used to manufacture commercial products without prior written approval from DVS Sciences Inc.

