

# LP 48.48 Control Line Fluid

## Purpose

The purpose of this reagent is to support library preparation with the LP 48.48 IFC (integrated fluidic circuit).

## Contents

Contents depend on product:

Part Number	Product Name	Quantity	Unit	Uses/Unit	Storage
101-2345	LP 48.48 Control Line Fluid—5 Pack	10 syringes	1 box	5 IFCs	Room temperature
101-2344	LP 48.48 Control Line Fluid—10 Pack	20 syringes	1 box	10 IFCs	

## Handling

Before using reagents, see the Advanta™ NGS Library Preparation with Juno™ Protocol (PN 101-7878) or the Advanta NGS Library Preparation with Access Array™ Protocol (PN 101-7885).

## Safety

Use standard laboratory safety protocols. Read and understand the safety data sheets (SDSs) before handling chemicals. To obtain SDSs, go to [fluidigm.com/sds](http://fluidigm.com/sds) and search for the SDS using either the product name or the part number.

**For technical support visit [fluidigm.com/support](http://fluidigm.com/support).**

North America +1 650 266 6100 | Toll-free (US/CAN): 866 358 4354 | [techsupport@fluidigm.com](mailto:techsupport@fluidigm.com) Latin America +1 650 266 6100 | [techsupportlatam@fluidigm.com](mailto:techsupportlatam@fluidigm.com)  
Europe/Middle East/Africa/Russia +44 1223 859941 | [techsupporteurope@fluidigm.com](mailto:techsupporteurope@fluidigm.com) China (excluding Hong Kong) +86 21 3255 8368 | [techsupportchina@fluidigm.com](mailto:techsupportchina@fluidigm.com)  
Japan +81 3 3662 2150 | [techsupportjapan@fluidigm.com](mailto:techsupportjapan@fluidigm.com) All other Asian countries/India/Australia +1 650 266 6100 | [techsupportasia@fluidigm.com](mailto:techsupportasia@fluidigm.com)

**For Research Use Only. Not for use in diagnostic procedures.**

Information in this publication is subject to change without notice. **Patent and license information:** [fluidigm.com/legalnotices](http://fluidigm.com/legalnotices). **Trademarks:** Fluidigm, the Fluidigm logo, Access Array, Advanta and Juno are trademarks and/or registered trademarks of Fluidigm Corporation in the United States and/or other countries.  
© 2018 Fluidigm Corporation. All rights reserved. 05/2018