


O-Ring, Nebulizer Adaptor Internal

Catalog#: 101817
Package Size: 5 units
Instrument: CyTOF and CyTOF 2

Technical Information

Description:

The Nebulizer Adaptor Internal O-ring is made of silicone (size number 2-009). The Nebulizer Adaptor connects the nebulizer to the spray chamber on the sample introduction system. It has two internal O-rings for sealing the contact of the adaptor and the nebulizer, and two external O-rings for sealing the contact of the adaptor and the spray chamber.



Contact Information:

Sales: sales@DVSsciences.com | Support: Support.NorthAmerica@DVSsciences.com or Support.Europe@DVSsciences.com
www.DVSsciences.com | For assistance by phone: +1-855-DVS-CYTO

This product is intended for Research Use Only. This product is not intended for therapeutic or diagnostic purposes in humans or animals.

Vendor info & description quote as needed.

List supplier of raw material antibody here only if targets/clones are unique, licensed, and/or indicated as required in our agreement.
© 2011 DVS Sciences Inc. The trademarks mentioned herein are the property of DVS Sciences Inc. or their respective owners.

DVS Sciences Inc. products are warranted to meet stated product specifications and to conform to label descriptions when used and stored properly. Unless otherwise stated, this warranty is limited to 6 months from date of sales for products used, handled and stored according to DVS Sciences, Inc. instructions. DVS Sciences, Inc. sole liability for the product is limited to replacement of the product or refund of the purchase price. Product may not be resold, modified for resale or used to manufacture commercial products without prior written approval from DVS Sciences Inc.

