

Torch Body

Catalog#: 101792
Package Size: 1 unit
Instrument: CyTOF and CyTOF 2

Technical Information

Description:

The CyTOF® Torch body is the part of the sample introduction system within which the plasma is ignited. The Torch body is a single-piece fused assembly of two concentric quartz tubes. When assembled, the Torch Body mounts onto the Torch Holder and is connected to the plasma gas and auxiliary gas lines. The Injector sits inside the torch body.



Contact Information:

Sales: sales@DVSSciences.com | Support: Support.NorthAmerica@DVSSciences.com or Support.Europe@DVSSciences.com
www.DVSSciences.com | For assistance by phone: +1-855-DVS-CYTO

This product is intended for Research Use Only. This product is not intended for therapeutic or diagnostic purposes in humans or animals.

Vendor info & description quote as needed.

List supplier of raw material antibody here only if targets/clones are unique, licensed, and/or indicated as required in our agreement.

© 2011 DVS Sciences Inc. The trademarks mentioned herein are the property of DVS Sciences Inc. or their respective owners.

DVS Sciences Inc. products are warranted to meet stated product specifications and to conform to label descriptions when used and stored properly. Unless otherwise stated, this warranty is limited to 6 months from date of sales for products used, handled and stored according to DVS Sciences, Inc. instructions. DVS Sciences, Inc. sole liability for the product is limited to replacement of the product or refund of the purchase price. Product may not be resold, modified for resale or used to manufacture commercial products without prior written approval from DVS Sciences Inc.

